

Gas Linings & Fume Hoods for AAS/ICP/MS/GC



We undertake laboratory gas handling and control systems for various inert, high purity, flammable and non flammable gases. To meet the requirements different industry application, we supply and install an extensive range of Gas Purification Panel. Our range includes Gas Purification and Control Panel, Cylinder Regulator Tubing and casing. We provide Gas Purification and Control Panel for 10/kg/cm inlet pressure & 7 kg/cm outlet pressure. Made using quality tested raw material. We undertake laboratory gas handling and control systems for various inert, high purity, flammable and non flammable gases. To meet the requirements different industry application, we supply and install an extensive range of Gas Purification Panel. Our range includes Gas Purification and Control Panel, Cylinder Regulator Tubing and Casing. We provide Gas Purification and Control Panel for 10/kg/cm

inlet pressure & 7 kg/cm outlet pressure. Made using quality tested raw material. Our Cylinder Regulator are of two types one in brass chrome plated body with SS diaphragm and other in SS body with SS diaphragm. The regulator comprises inlet pressure gauge displays the total pressure of cylinder and outlet pressure gauge shows the regulated pressure. We supply Tubing Casing, which is available in 1/2", 1/4", 1/8", 1/16" sizes that comes in copper, S.S. Teflon, nylon material. The size of the tubing casing depends upon the required flow and pressure. Based on cutting edge technology, the casing capping is used to cover the tubing and to maintain the aesthetic look. We offer three different types of casing capping, which are available in M.S. powder coated, aluminum powder coated and PVC material form.