

Laboratory Furnitures

All types of Laboratory Furniture's like Island Tables, Instrument Tables, Wall Tables, Anti Vibration Tables, Wall Cupboards, Chemical Storage Cabinets, Reagent Racks, SS Stools. Custom made lab furniture consists of under-bench modular with inner storage shelves, lockable doors, granite working tops, foot rest, reagent racks with power points, water sinks and inlets & outlets etc.

Material of Construction: All Lab Furniture's are available in the following materials of construction depending upon customer's choice.

- Wooden laminated finish
- The Tables are provided with Top Drawer and Bottom Cupboard with shelf and shutter mounted on imported SS spring loaded hinges. The handles are of SS material. The Drawers and Shutters are provided with locks.
- PP / SS / Porcelain Sinks with Swan Neck Taps and Electrical Power Modules can be provided for all tables as options.
- Table size and colour as per the customer requirement. Custom tables to meet customers / work area requirement can be offered.

