

Muffle Furnace



Quesst International

Light weight model with ceramic insulation grooved type outer body made of thick mild steel sheet. Inside chamber is made of silluminate grooved refractory for supporting the heating elments on four sides. Heating elements are made of imported kanthal wire and easily replaceable. Two working temperature of 960 & 1000 can be provided as per requirment.complete with digital temerature controller cum indicator.indicator, thermocouple, copper thermal fuse are also provided. Works on 200/220 volts ac.

Temperature Range	Inner Dimension in mm			Rating(kw)
Upto 930 C Upto 1000 C	100	100	230	1.75
	130	130	250	2.5
	150	150	300	3
	100	100	230	2
	130	130	250	3
	150	150	300	3.75