

Water Bath-Digital Temperature Controlled

This is a regular double walled water bath. The inner tank is made of ss and the outer is Made of mild steel powder coat painted. Stirrer of 1/35 hp motor keeps the temperature Uniform. The temperature is controlled By digital temperature controller cum indicator. An ss tray is provided at the bottom for holding the samples.

SIZE: 400MM L X 300MM B X 100MM H TEMP.RANGE: 5DEG. C ABOVE AMBINT TO100 DEG. C ACCURACY: ± 1 DEG.C
--

ACCURACY: ± 1 DEG.C

(ANY OTHER SIZE CAN ALSO BE MANUFACTURED UPON REQUEST)

Quesst International