

## Water Bath Incubator Shaker-(Metabolic Shaking Incubator)



Temperature range from room temperature to 100 °c.digital temperature control with an accuracy of plus minus 0.5 °c. double walled inside made of stainless steel 304 qlty and outside mild steel sheet powder coat painted. Oscillating tray is riding on ball rollers. Shaking speed range 40 to 140 r.p.m.built in water drain and digital rpm indicator to work on 220/230 volts a.c. shaking tray size : inside (stainless steel 304 qlty)

SHAKING TRAY SIZE : INSIDE (STAINLESS STEEL 304 QLTY)
---

Isb-1 275 x 275 x 75 mm
-------------------------

Isb-2 275 x 275 x150 mm
-------------------------

Isb-3 405 x 300 x 150 mm
--------------------------

SHAKING RACKS MADE OF STAINLESS STEEL FOR TYPE ISB 1 & ISB 2
--

To hold 16 nos. 25 ml conical flask
-------------------------------------

To hold 12 nos. 50 ml conical flask
-------------------------------------

To hold 9 nos. 100 ml conical flask
-------------------------------------

To hold 5 nos. 250 ml conical flask
-------------------------------------

To hold 4 nos. 500 ml conical flask
-------------------------------------

To hold 2 nos.1000 ml conical flask
-------------------------------------

To hold 25 nos. 25 mm dia test tube
-------------------------------------