

Water Bath-Thermostatic Controlled(Double Walled)



Temperature range from ambient to 100 deg.c. electrically heated model thermostatic control with an accuracy of plus or minus 2 deg.c. double walled inside stainless steel and outside mild steel powder coat painted Top of the bath is provided with a cover made of aluminium powder coat painted works on 220 volts ac Inner water tank is made of 304 ss sheet and argon arc welded and hence leak proof for a minimum period of 5 years unlike other tanks which are of light gauge and lead soldered. After a few operations leak may occur spoiling glass wool and thermostat. Concentric rings are of 304 ss quality .

| Description | Size in mm | Rating (kw) |
|---------------------|-----------------|-------------|
| 6Holes of 75mm dia | 300 x 250 x 100 | 1 |
| 12Holes of 75mm dia | 355 x 355 x 405 | 1.5 |

Quesst International